



An ISO 9001 - 2008 Certified Company

ADI DUO PLUS

Semi Auto Chemistry Analyzer With Coagulation



SALIENT FEATURES

- 7 Inch LCD Touch Screen Display
- Advance Optics Filter System
- 7 Filters 340, 405, 450, 505, 546, 578, 630
- Mode of Testing : End Point, Fixed Time, Kinetic
- Open Reagent System
- Low Carry Over
- Reaction Time Display
- Quality Control for all Test
- Add or Delete More Than 300 Programs Easily
- User Friendly Interface
- Universal Power Supply

TEST PARAMETERS

ADA, D-HDL, D-LDL, ASO, CRP, RF, Microalbumin, GLU, TG, Chol, Urea, Uric Acid, SGOT, SGPT, ALP, ALB, Creatinine, Bilirubin T/D, Calcium, LDH Sodium, Potassium, Magnesium, Iron, Phosphorus, Amylase, GAMMA GT, TP, G6PDH, Ammonia, CK-NAC, CK-MB.

All Clinical Chemistry Reagents, and Turbidometry Kits Can be done.

ADI DUO PLUS

Semi Auto Chemistry Analyzer With Coagulation



Technical Specifications


- Analysis Method: 1 Point End, 2 Point End (Sample Blank) Kinetics, Fixed-Time.
- Optics: 7 Filters, 340, 405, 450, 505, 546, 578, 630.
- Lamp House: Long-life Tungsten Halogen Lamp (6V 10W) Auto-Sleep Function
- Display 7.0" TFT LCD with Touch Screen, 800*480 Pixels
- Printer: Build-in Thermal Printer, 57mm Printer Paper
- Temperature: 25, 30, 37°C and Room Temperature.
- Flow Cell Capacity & Optical Path: 32ul, 10mm
- Resolution: 0.0001Abs
- Photometric Range: 0.0000-3.0000 Abs
- CPU: ARM Cortex-A8, 720MHz
- Memory: DDR2 SDRAM: 1G bit
- Flash: 256 MB
- Over 300 Programmers, 200000 Results can be Stored
- Interface: 4 USB Host; 1 USB Slave; 1 Ethernet Port
- USB Ports (Host) for Keyboard, Mouse, Printer, U Disk Bar Code Reader etc.
- USB Ports (Slave) for PC Connection (Optional) (Support PC Software Control Function)
- Environment: Temperature: 10°C~37°C Humidity: =85%
- Power Requirement: Wide Power Supply AC100 240V, 50/60Hz
- Dimension: 400 x 300 x 160mm (LXWXH)
- Weight: 7.5kg



MAIN SCREEN



TEST SCREEN



PARAMETER SETUP



For Further Details Contact Us :



ADI DIAGNOSTICS PVT. LTD.

Head Office: 11-5-422/1, Flat No. 301, Kiran Apartments, Red Hills, Lakdi-ka-pool, Hyderabad - 500 004, Telangana State, India.
Mob: 9848018523, 9848195080,
E - Mail: info@adidiagnostics.com

Factory: Plot No. L - 31, Green Industrial Park, TSIC, Polepally SEZ, Jadcherla - 509 301, Mahboobnagar Dist, Telangana State, India.
Website: www.adidiagnostics.com

Distributor

Specifications subject to change without prior notice